

Mebiol®

Thermoreversible Hydrogel

Catalog #5180-10ML

Full Product Description

Mebiol® Thermoreversible Hydrogel is a copolymer of poly(N-isopropylacrylamide) and poly(ethylene glycol) (PNIPAAm-PEG). The defining feature of Mebiol® is its temperature reversible sol-gel transition. Chill Mebiol® to turn it into a sol (meaning it handles like a liquid). Heat Mebiol® to turn it into a solid hydrogel. The unique thermoreversible feature allows for extremely easy cell handling. Cultures are seeded into cooled Mebiol® and recovered conveniently by cooling the culture vessel and centrifugation. In the gel state, the highly lipophilic and network structure of the Mebiol® presents an efficient niche for cell proliferation, cell communication, gas and mass exchange, and protection of cells and tissue from shear forces.

Mebiol® is easy to use, non-toxic, 100% synthetic and has high transparency. Published applications include stem cell and pluripotent stem cell culture, expansion, differentiation, spheroid culture, cell implantation, organ and tissue regeneration, drug delivery and more. Mebiol® has been commercialized as a cell/tissue culture reagent for ES cells, chondrocytes and cancer cells.

Parameter, Testing, and Method	Mebiol® Thermoreversible Hydrogel Catalog #5180-10ML
Package Size	10 mL
Form	Lyophilized
Storage Temperature	Room Temperature (25° C)
Fluid Temperature	0°-15° C
Gel Temperature	>25° C
Opacity	High Transparency
Sterilization method	Ethylene Oxide

Preparation and Usage

Cell culture in Mebiol® flask:

1. Open package in sterile environment and add 10mL of culture media directly to Mebiol® flask.
2. Close flask and lay in refrigerator (2-10 °C) for 3 hours to allow Mebiol® to absorb media.
3. Periodically shake flask gently to help dissolve the Mebiol® into the media. Make take up to a day for complete dissolution (maintain at cold temperature).
4. After dissolving, allow to sit in refrigerator until bubbles are eliminated.
5. Add cells/tissues into Mebiol® at low temperature (2-10°C) and warm to 37°C to culture cells/tissues in a 3D hydrogel environment.

6. To recover the cells/tissues, cool the Mebiol® back down and dilute with 30-40 mL of cold saline or medium. This will take it out of gel phase. Suspended cells can be easily recovered via centrifugation.

Cell culture in 24-well plate:

1. Follow steps 1-4 above.
2. Before adding cells/tissues, remove desired amount of Mebiol® and put into chilled centrifuge tube. Put remaining Mebiol® into refrigerator (still in original flask).
3. Add approximately 10^5 cell/mL to chilled centrifuge tube and mix gently on ice.
4. Warm up 24-well plate and separate overlaying culture media to 37°C.
5. Add 200 μ L of cold Mebiol® suspension to the center of each well. Make sure Mebiol® does not cover entire bottom of the well to make it easier to change overlaid medium in later step.
6. Overlay 400 μ L of warmed culture media containing phenol red on Mebiol® suspension. Culture at 37°C.

7. To recover cells, refrigerate at 2-8°C and dilute with 400 μ L of cold saline or medium. This will take Mebiol® out of gel phase. Suspended cells can be easily recovered via centrifugation.

Do not re-sterilize to avoid deterioration. Add media promptly after opening, store in refrigeration and use within 1 month.

Application Examples

Cell Proliferation – Cancer tissue grown in 2D or collagen can be overrun by fibroblasts. In Mebiol®, fibroblast growth is suppressed allowing proliferation of cancer cells.

Formulation of Spheroid – MEC (Muco-Epidermoid Carcinoma) was cultured in Mebiol® for 3 weeks and a spheroid formed.

Cell Preservation – Allows you to cultivate human cells long-term while maintaining their structure and integrity.

Disclaimer

This product is for R&D use only and is not intended for human or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.